

## Baby Poser Instructions

### Thank you for purchasing a Lastolite By Manfrotto Baby Poser

This innovative accessory has been designed especially for baby photography and consists of three main sections (see diagrams). These can be arranged in a number of ways to facilitate easy posing of even the smallest baby.

**Contents:** Main cushion; End cushion/neck support; Centre cushion/seat; velour throwover background (doubles as a carry bag).

#### 1. Baby poser in the flat position


Used for posing babies either on their backs or their tummies.

All three sections are velcroed together in the flat position, using the velcro joining tapes.

The adjustment straps can be velcroed together under the unit, tightly to loosely to assist moulding or flattening of units.


The central section may be removed to accommodate a bulky nappy and reduce its prominence in the photograph.

A small cushion or pad may be placed under the unit at the head end to increase the slope.



#### 2. Poser sections positioned in chair formation

The end cushion is placed on top of the main cushion to form a chair. The center cushion is fitted into the central gap for added support of left out completely for accommodation of baby's bottom.



#### 3. Care of your baby poser

The velour drape cover may be dry-cleaned if it becomes necessary.